

BIO-BUFFER

ST08006-125/500

DESCRIPTION

1X filtered PBS solution pH 7.4 ready-to-use. This solution does not contain calcium or magnesium. It has similar physicochemical conditions as the human body, such as pH, osmolarity and ion concentrations.

RECOMMENDED USES

To dilute, wash or resuspend biological samples of cell and tissue cultures, or antioxidant standards. Compatible with e-BQC lab and e-BQC lab CEAC devices.

COMPATIBLE STANDARDS

Product	Reference
Ascorbic acid standard	ST08001
Gallic acid standard	ST08002
Trolox standard*	ST08003

*When diluting Trolox, check the pH. If lowered, use NaOH to keep it at 7.4 to facilitate solubility.

ADDITIONAL INFORMATION

Designed for research use only, it is not approved for human or animal use, or clinical diagnosis.

Storage: Room temperature in a dry place

Best before: 1 year from the date of delivery

